



IWeb Release Notes

v5.20.0



Support Services

For general support on this product, contact your system administrator or help desk. For up-to-date documentation, visit the STC Documentation Portal at <https://documentation.stchome.com/>.

Connect with Us on Social Media



Copyrights and Trademarks

© 2019 by Scientific Technologies Corporation (STC). All rights reserved.

This documentation describes the following: IWeb (v5.20.0) release notes

No part of this publication may be altered, reproduced, transmitted, stored in a retrieval system, or translated into any human or computer language, in any form or by any means, electronic, mechanical, magnetic, optical, or otherwise, without the prior written permission of the copyright owner.

This document, along with any embedded media, is provided for informational purposes only. Scientific Technologies Corporation (STC) provides this material "as is" and does not assert that this document is error free. The information contained herein may not necessarily accurately represent the current or planned functions of the application, and may be subject to significant and frequent modification. Scientific Technologies Corporation (STC) may change this document, any embedded content, and the product described herein at any time. Any changes will be incorporated in new versions of this document.

Scientific Technologies Corporation (STC) recognizes the rights of the holders of all trademarks used in its publications.

This document may provide hyperlinks to third-party websites or access to third-party content. Links and access to third-party sites are provided for your convenience only. Scientific Technologies Corporation does not control, endorse, or guarantee third-party content and is not responsible for any content, associated links, resources or services associated with a third-party site. Scientific Technologies Corporation shall also not be liable for any loss or damage associated with your use of any third-party content.

Table of Contents

Support Services.....	2
Connect with Us on Social Media	2
Copyrights and Trademarks	2
Table of Contents.....	3
New Features	4
Fixed Bug List	5
Known Issues.....	5
Fixed Bug Details.....	5
Export to Excel.....	9
Product Documentation.....	10
Database Release Notes	10

New Features/Improvements

The following are new features, functionalities or enhancements in this version of IWeb.

Key	Summary	Description
CONSORT-4911	State field is required on Facility maintenance	<p>Made changes to have 'State' as a required field on Facility Maintenance page.</p> <p>Steps to Test:</p> <ol style="list-style-type: none"> 1. Login to IWeb as a Registry Client user. 2. Select an existing Organization. 3. Go to Facilities > Search Add. 4. Click on Search button. 5. Click on 'Add' button to create a new Facility. 6. Notice the State drop-down is displayed with the current logged in state code as default value. 7. Notice that there is no Blank option in the state drop-down.
CONSORT-4893	As an administrator, I want to be able to have VFC status as required, so that users are required to enter the information	<p>Added new functionality where a VFC Status is required when adding a facility.</p> <p>Steps to Test:</p> <ol style="list-style-type: none"> 1. Login as a Registry Client. 2. Select an organization. 3. Go to Facility > Search / Add. 4. Start to add a new facility. 5. Notice the VFC Status field appears in red, indicating it is a required field. 6. Leave the VFC Status field blank, complete the required fields, and save. 7. Notice a red message is displayed and the facility is not saved: VFC Status must be entered. 8. Enter a VFC status and save. 9. Verify the facility is saved without error. 10. Click on the Edit button, clear the VFC Status and save. 11. Notice the edit is prevented, and a VFC status must be selected to save.

Fixed Bug List

The following bugs were fixed in this version. For detailed information, see the [Fixed Bug Details](#) section below.

Key	Summary
CONSORT-4908	Remove serology module from the AZ build
CONSORT-4909	SSO - Properties page edits settings that aren't displayed, blanking them out.
CONSORT-4910	Separate Bad Merge: Guardian details are not displayed correctly
CONSORT-4912	WA - Health Promotions Message Change
CONSORT-4920	SSO - Default password settings on IWeb is not working correctly when connected
CONSORT-4921	Performance issue on Vaccine Administered Report
(no steps to test)	Security work – Priority High – SQL Injections

Known Issues

There are no known issues with this release.

Fixed Bug Details

The following table lists the detailed information about each of the bugs fixed in this version. To export the testing steps to Excel, see the [Export](#) section below.

Key	Summary	Affects Client	Description
CONSORT-4908	Remove serology module from the AZ build	AZ	<p>Steps to Test:</p> <ol style="list-style-type: none"> 1. Start with a 5.19.0 or earlier version of IWeb (an earlier version that had Serology turned on.) 2. Go to Administration > Properties > Registry Settings. 3. Notice the Enable Serology property is on the page. 4. Go to Administration > Config Serology. 5. Notice the page opens and allows serology codes to be added. 6. Select an organization. 7. Go to Patient> Search / Add and select a patient. 8. In the left side navigation bar, go to Serology.

Key	Summary	Affects Client	Description
			<ol style="list-style-type: none"> 9. Notice the page opens and allows a serology record to be opened. 10. Upgrade to IWeb v5.20.0+. 11. Go to Administration > Properties > Registry Settings. 12. Notice the Enable Serology property not on the page. 13. Go to Administration > Config Serology. 14. Notice the menu option is not on the page. 15. Select an organization. 16. Go to Patient> Search / Add and select a patient. 17. In the left side navigation bar, notice Serology is not an option.
CONSORT-4909	SSO - Properties page edits settings that aren't displayed, blanking them out.	All	Fixed issue where, in a SSO environment, checkboxes could be disabled accidentally on the properties page. Setup: <ul style="list-style-type: none"> • System has SSO and is on IWeb 5.19.2.0+. <p>Steps to Test:</p> <ol style="list-style-type: none"> 18. Login as a Registry Client user. 19. Go to Administration > Properties. 20. Enable a checkbox. 21. Save. 22. Verify the url ends with editProperties.do 23. Click on the address bar. 24. Press Enter. 25. Go back to the properties page. 26. Notice the enabled checkboxes are still enabled.
CONSORT-4910	Separate Bad Merge: Guardian details are not displayed correctly	All	Fixed the issue where Guardian details are not displayed correct after Separate Bad Merge, when Multi-rows property is turned off for Family and Contact Information. Setup: <ul style="list-style-type: none"> • In Administration > Properties > Registry Settings > Enable users to report duplicates is enabled. • Have Multi-rows property un-checked under Family and Contact Information.

Key	Summary	Affects Client	Description
			<p>Steps to Test:</p> <ol style="list-style-type: none"> 1. Login as a Registry Client user. 2. Select an Organization/Facility (Organization A/ Facility A). 3. Go to Patient Search/Add. 4. Create a Patient with guardian added (Patient A). 5. Go to Patient > Search/Add. 6. Create another Patient with same date of birth and with guardian details added (Patient B). 7. Go to Patient > Search / Add and Search by date of birth of above patients. 8. Notice that Patient A and Patient B are returned in Search Results. 9. Click on the Report Duplicates button and check the check-box corresponding to both Patient A and Patient B. 10. Select a Reason for deduplication and Click 'Report Duplicates' button. 11. Select Patient B as master and Click on 'Merge'. 12. Go to Administration > Manual Deduplication. 13. Select the 'User Flagged' radio button and click on 'Continue' button. 14. Scroll through the records until Patient A and Patient B are displayed. 15. Click on the Merge button. 16. Go to Administration > Separate Bad Merges. 17. Enter the SIIS ID of Patient B and Click on 'Get Reserve Mapping Records' button. 18. Select Patient B as master and Select reserve records to separate. 19. Click on 'Separate Reserve Records' button. 20. Notice that message is displayed as Records Separated. 21. Go to Patient > Search/Add. 22. Search for the above 2 Patients. 23. Select and open both the Patients to view Patient demographic page. 24. Notice that Patient A and B has the correct guardian details displayed.

Key	Summary	Affects Client	Description
CONSORT-4912	Health Promotions Message Change	WA	<p>Modified Health Promotions so that a single confirmation process is used to set a patient's valid email address to Yes. Also, modified the thank you message.</p> <p>Setup:</p> <ul style="list-style-type: none"> Registry Client user has Health Promotions Administrator permission. Test system interfaces to Health Promotions. <p>Steps to Test:</p> <ol style="list-style-type: none"> Login as a Registry Client user. Select an Organization. Go to Patient > Search Add. Select or a patient who is six years or less old. Patient should have first name, last name, date of birth, guardian first name, guardian last name, and an address with a WA zip code. Patient should not have an email address to start. It is a good idea to note or screen capture the patient's information. Click on the Health Promotion Info button. Click on the Update button. Notice Valid Email Address is set to N. Access Health Promotions and start to enroll the patient in email. Enter an email address for the patient and complete the required fields so they are an exact match of the patient information. Submit. Verify the confirmation message reads correctly in the English and Spanish version.
CONSORT-4920	SSO - Default password settings on IWeb is not working correctly when connected	All	<p>Made change so that an error message is displayed if Default Password cannot be re-used.</p> <p>Setup:</p> <ul style="list-style-type: none"> Have the Password Rules set in Keycloak > Authentication > Password Policy page and IWeb> Administration > Password Rules in sync. <p>Steps to Test:</p> <ol style="list-style-type: none"> Login to IWeb as a Registry Client user with System Admin Permissions.

Key	Summary	Affects Client	Description
			<ol style="list-style-type: none"> 2. Go to Administration > Search/Add user. 3. Search and Click Add button. 4. Check the 'Use the default system password "XXXXXXXXXX!" and Expire Now check-box and Create a user by entering all the Required details. 5. Logout. 6. Login to the system with the credentials of user created in Step 4. 7. Enter the New Password on Change Password screen and Submit. 8. Logout. 9. Repeat Steps 1 to 2. 10. Search for and Select the user created in Step 4. 11. Click Update. 12. Check the 'Use the default system password "XXXXXXXXXX!" check-box to reuse the default password. 13. Click on Update. 14. Notice that error message is displayed as: Cannot re-use a password that has been used in the past XX passwords. Please try another password. (XX represents the number of passwords that cannot be reused as per the system and SSO password configuration).
CONSORT-4921	Performance issue on Vaccine Administered Report	All	<p>Made changes to improve the performance when Vaccine Administered Report is run.</p> <p>Steps to Test:</p> <ol style="list-style-type: none"> 1. Login to IWeb as a Registry Client user. 2. Go to Reports > Report Module > Vaccine Administered Report. 3. Click on Create Report with default parameters (i.e Do not Limit). 4. Notice that a new window is opened with report results.

Export to Excel

Follow these steps to export the testing steps to an Excel spreadsheet:

1. Click this link and log in to Jira if required:
<https://stchome.atlassian.net/issues/?filter=52665>.

2. Click the **Change View** icon and select **List View**.
3. Click the **Export** icon (it looks like a download icon) and export as needed.

Product Documentation

Product documentation is located on the STC Documentation Portal:
<https://documentation.stchome.com/>.

The following documents are available on the Documentation Portal for this release of IWeb:

- IWeb (v. July 2018) User Guide (no changes for this version)
- IWeb (v. March 2018) Quick Reference Guides (no changes for this version)

Database Release Notes

The following briefly addresses changes in the database for users with read-only access:

- Added VFC_STATUS_CODE to H33_IRMS Table.
- Added pre booking DB changes
- Fixed issue to show the right guardian details on separating a bad merge
- Set enableSerology on H33_KEYS to N for AZ
- [CONSORT-4907](#) - Remove SEROLOGY_MANAGEMENT_ACCESS, SEROLOGY_VIEW_ACCESS permission from all AZ users
- [CONSORT-4918](#) - Modified vaccination reserve audit trigger to only set Deletion on actual removal of data
- Fixed vaccine audit trigger to populate values in the audit columns on delete