





IWeb Release Notes



Support Services

For general support on this product, contact your system administrator or help desk. For up-to-date documentation, visit the STC Documentation Portal at https://documentation.stchome.com/.

Connect with Us on Social Media



Copyrights and Trademarks

© 2019 by Scientific Technologies Corporation (STC). All rights reserved.

This documentation describes the following: IWeb (v5.30.0) release notes

No part of this publication may be altered, reproduced, transmitted, stored in a retrieval system, or translated into any human or computer language, in any form or by any means, electronic, mechanical, magnetic, optical, or otherwise, without the prior written permission of the copyright owner.

This document, along with any embedded media, is provided for informational purposes only. Scientific Technologies Corporation (STC) provides this material "as is" and does not assert that this document is error free. The information contained herein may not necessarily accurately represent the current or planned functions of the application, and may be subject to significant and frequent modification. Scientific Technologies Corporation (STC) may change this document, any embedded content, and the product described herein at any time. Any changes will be incorporated in new versions of this document.

Scientific Technologies Corporation (STC) recognizes the rights of the holders of all trademarks used in its publications.

This document may provide hyperlinks to third-party websites or access to third-party content. Links and access to third-party sites are provided for your convenience only. Scientific Technologies Corporation does not control, endorse, or guarantee third-party content and is not responsible for any content, associated links, resources or services associated with a third-party site. Scientific Technologies Corporation shall also not be liable for any loss or damage associated with your use of any third-party content.

Table of Contents

Support Services	2
Connect with Us on Social Media	
Copyrights and Trademarks	2
Table of Contents	3
New Features/Improvements	4
Fixed Bug List	4
Product Documentation	5
Database Release Notes	5

New Features/Improvements

The following are new features, functionalities or enhancements in this version of IWeb.

Кеу	Summary	Affects Client	Description
	Patients	All	Increased Patient's Last Name field character limit to 100 characters.
	DTT	All	Enhancement to allow designated providers to import inventory using the DTT. In order for a provider to have the ability to import inventory, the property 'Import Inventory Data From IIS' needs to be enabled on the Organization or Facility Maintenance screen.
	Facility Merge	All	 Facility Merge Logs will display on the Organization / Facility Merges page. The new section will include the below details: Merge Date Site Eliminated Destination Site Merged By Start Date End Date Merge Error 'Refresh View' button will refresh the table and display all the Facility Merges that are done and currently in progress.
HDSD-1594	Required Fields	All	 Added the ability to mark the following fields as required when adding and editing Patient information. Patient Phone Number Patient Nationality
<u>HDSD-1542</u> and <u>HDSD-</u> <u>1395</u>	Provider Agreement	TN	Updated the TN Provider Agreement PDF to display the new CDC language.

Fixed Bug List

The following table lists the detailed information about each of the bugs fixed in this version.

Кеу	Summary	Affects Client	Description & Impact
HDSD-1279 HDSD-138 HDSD-265	Approved Orders		Previously, order items would randomly be duplicated during the ordering process. Now, order item will not be duplicated because the Approve Orders button will be disabled after the first selection.
HDSD-1930	IMM Cert PDF	TN	Previously, the IMMCert PDF could be edited while opened

Кеу	Summary	Affects Client	Description & Impact
			in Chrome. This has been resolved and end-user edits will no longer be an option.
	DTT	ALL	The DTT will now accept a date format for the patient date of birth with a timestamp even when the date format is set without one in the profile. For example, when a profile is defined with a patient DOB as DD/MM/YYYY, the system will now accept data imported as DD/MM/YYYYHH:mm:ss and will strip off the timestamp in order to insert the data properly.

Product Documentation

Product documentation is located on the STC Documentation Portal: <u>https://documentation.stchome.com/</u>.

The following documents are available on the Documentation Portal for this release of IWeb:

- IWeb (v. July 2018) User Guide (no changes for this version)
- IWeb (v. March 2018) Quick Reference Guides (no changes for this version)

Database Release Notes

The following briefly addresses changes in the database for users with read-only access:

 Code changes to add new match set 125(PASSPORT+DOB) in H33_WEIGH, corresponding changes to evaluate and pass PASSPORT in H33_REPLACE_CHECK and H33_FIND_DUPLICATESUNO-2559. Increased column length of patient last name to 100.