



# IWeb Release Notes

v5.32.0



## Support Services

For general support on this product, contact your system administrator or help desk. For up-to-date documentation, visit the STC Documentation Portal at <https://documentation.stchome.com/>.

## Connect with Us on Social Media



## Copyrights and Trademarks

© 2019 by Scientific Technologies Corporation (STC). All rights reserved.

This documentation describes the following: IWeb (v5.32.0) release notes

No part of this publication may be altered, reproduced, transmitted, stored in a retrieval system, or translated into any human or computer language, in any form or by any means, electronic, mechanical, magnetic, optical, or otherwise, without the prior written permission of the copyright owner.

This document, along with any embedded media, is provided for informational purposes only. Scientific Technologies Corporation (STC) provides this material "as is" and does not assert that this document is error free. The information contained herein may not necessarily accurately represent the current or planned functions of the application, and may be subject to significant and frequent modification. Scientific Technologies Corporation (STC) may change this document, any embedded content, and the product described herein at any time. Any changes will be incorporated in new versions of this document.

Scientific Technologies Corporation (STC) recognizes the rights of the holders of all trademarks used in its publications.

This document may provide hyperlinks to third-party websites or access to third-party content. Links and access to third-party sites are provided for your convenience only. Scientific Technologies Corporation does not control, endorse, or guarantee third-party content and is not responsible for any content, associated links, resources or services associated with a third-party site. Scientific Technologies Corporation shall also not be liable for any loss or damage associated with your use of any third-party content.

# Table of Contents

<b>Support Services.....</b>	<b>2</b>
<b>Connect with Us on Social Media.....</b>	<b>2</b>
<b>Copyrights and Trademarks.....</b>	<b>2</b>
<b><i>Table of Contents .....</i></b>	<b>3</b>
<b><i>New Features/Improvements .....</i></b>	<b>4</b>
<b><i>Fixed Bug List .....</i></b>	<b>4</b>
<b><i>Product Documentation.....</i></b>	<b>5</b>
<b><i>Database Release Notes .....</i></b>	<b>5</b>

# New Features/Improvements

The following are new features, functionalities or enhancements in this version of IWeb.

Key	Summary	Affects Client	Description
	Security	All	Updated the code to remove security vulnerability in the following locations: <ul style="list-style-type: none"> <li>• Allocations</li> <li>• Daily Inventory</li> <li>• Child Care</li> <li>• Lot Numbers</li> <li>• Vaccination Data Quality Report</li> <li>• First Responder</li> <li>• Health Promotions (WA)</li> <li>• Patient Data Quality Report</li> <li>• Default Settings</li> <li>• Patient</li> </ul>
	Serology	All	Added the ability for Serology tests to be entered using HL7 through PHC Hub. These tests are deduplicated and are viewable in from the user interface. Detailed document to follow.
<a href="#">HDSD-1927</a>	Anonymous	LA	Enabled SIIS Opt Out and Add Anonymous Vaccinations for LINKS
<a href="#">HDSD-1374</a>	Certificate of Immunization	MS	The following changes have been made to the Certificate of Immunization for MS. <ul style="list-style-type: none"> <li>• Updated the validation for MMR on the 121 form.</li> <li>• Change the childcare certificate rules to show a temporary certificate until the child has met all requirements</li> </ul>
<a href="#">HDSD-1243</a>	Patient Record Report	All	Add a new property, 'Hide Physician Signature line from Patient Record Report'. When enabled the Physician Signature option will not display on the Patient Record Report.

## Fixed Bug List

The following table lists the detailed information about each of the bugs fixed in this version.

Key	Summary	Affects Client	Description & Impact
	Lot Numbers	All	Previously, a stack trace was displayed when trying to manually add a duplicate lot. Now, an error message is displayed stating the lot is not saved because the lot already exists.

Key	Summary	Affects Client	Description & Impact
	First Responder	All	Previously, the following sections of First Responder display error message or stack traces with navigating to the page. Now, the page will load without error. <ul style="list-style-type: none"> <li>• Reminder Recall</li> <li>• Patient Detail Report</li> <li>• Lot Numbers – Add Lot</li> </ul>
<a href="#">HDS-1831</a>	Quantity on Hand	All	Previously, the quantity on hand display on the Lot Number Selection module when administering a vaccine was incorrect. Now, the quantity on hand displays correctly.
<a href="#">HDS-1002</a>	PAIS	All	Previously, an incorrect code around organizational status was given. Now, the correct code is given.

## Product Documentation

Product documentation is located on the STC Documentation Portal: <https://documentation.stchome.com/>.

The following documents are available on the Documentation Portal for this release of IWeb:

- IWeb (v. July 2018) User Guide (no changes for this version)
- IWeb (v. March 2018) Quick Reference Guides (no changes for this version)

## Database Release Notes

The following briefly addresses changes in the database for users with read-only access:

- Code changes to add new match set 125(PASSPORT+DOB) in H33\_WEIGH, corresponding changes to evaluate and pass PASSPORT in H33\_REPLACE\_CHECK and H33\_FIND\_DUPLICATESUNO-3538 Remove all references to h33\_serology\_master
- Data migration for client's who currently use the serology module
- Created new table H33\_SEROLOGY\_INVALID\_DATA to hold invalid serology data
- Updated link serology to populate ASIIS\_PAT\_ID in serology reserve
- Created new lookup INVALID\_SEROLOGY\_REASON (for INVALID\_REASON in H33\_SEROLOGY\_INVALID\_DATA)
- Modified serology dedup to move invalid data to H33\_SEROLOGY\_INVALID\_DATA
- Added code in dedup to capitalize IRMS\_PAT\_ID (to match patient reserve) and other relevant fields in H33\_SEROLOGY\_PRE\_RESERVE table - needed to populate ASIIS\_PAT\_ID
- Change SPECIAL\_CONSIDERATION contra\_type to CONTRAINDICATION
- Correct text NULL in h33\_serology\_test\_result.result\_loinc\_code to NULL value
- PAIS - Made changes so Patient Reserve being saved correctly for Inactive Code

- Added mapping\_key column to h33\_pre\_booking\_request table and changed flu\_year to be nullable
- Removed additional code added to dedup