



# ImmuCast

## Release Notes

v5.23.0



## Support Services

For general support on this product, contact your system administrator or help desk. For up-to-date documentation, visit the STC Documentation Portal at <https://documentation.stchome.com/>.

## Connect with Us on Social Media



## Copyrights and Trademarks

© 2019 by Scientific Technologies Corporation (STC). All rights reserved.

This documentation describes the following: ImmuCast 5.23.0 (and IWeb Forecaster) release notes

No part of this publication may be altered, reproduced, transmitted, stored in a retrieval system, or translated into any human or computer language, in any form or by any means, electronic, mechanical, magnetic, optical, or otherwise, without the prior written permission of the copyright owner.

This document, along with any embedded media, is provided for informational purposes only. Scientific Technologies Corporation (STC) provides this material "as is" and does not assert that this document is error free. The information contained herein may not necessarily accurately represent the current or planned functions of the application and may be subject to significant and frequent modification. Scientific Technologies Corporation (STC) may change this document, any embedded content, and the product described herein at any time. Any changes will be incorporated in new versions of this document.

Scientific Technologies Corporation (STC) recognizes the rights of the holders of all trademarks used in its publications.

This document may provide hyperlinks to third-party websites or access to third-party content. Links and access to third-party sites are provided for your convenience only. Scientific Technologies Corporation does not control, endorse, or guarantee third-party content and is not responsible for any content, associated links, resources or services associated with a third-party site. Scientific Technologies Corporation shall also not be liable for any loss or damage associated with your use of any third-party content. (20181212)

# Table of Contents

Introduction.....	1
Apply the Release.....	1
For ImmuCast (Stand-Alone Forecaster).....	1
For IWeb Only.....	1
Patient Groups Affected by this Release (v5.23.0) .....	1
New Features/Improvements .....	2
Bug Fix Summary .....	2
Product Documentation .....	3

## Introduction

This release applies new CVX vaccine codes only. There are no CPT code updates. The new vaccine codes are not included in forecasting logic.

## Apply the Release

Apply the release by executing either the included `forecast.bat` or `forecast.sh` file. It is not necessary to verify which forecast releases have been previously applied; however, to ensure that the system is up to date with new CVX codes, prior releases should have already been applied. Log files, which are created in the same folder from where the release is executed, can be reviewed for errors. To determine the current version of ImmuCast, execute this statement from SQLPlus:

```
SQL> select max(version) from h33_forecast_version where insert_stamp =  
(select max(insert_stamp) from h33_forecast_version);
```

If the version number returned is not 5.22.3, you may download and apply the previous releases prior to applying this release. Log files are created in the folder from where the release is executed and can be reviewed for errors.

### For ImmuCast (Stand-Alone Forecaster)

After applying the release, restart Tomcat to enable and cache any new vaccine codes into memory.

### For IWeb Only

Since this release does not include forecasting changes, no reforecasting of patients is needed.

## Patient Groups Affected by this Release (v5.23.0)

- None

## New Features/Improvements

The following are new features, functionalities or enhancements in this version of ImmuCast.

Ticket			
Description			
	CVX Code	CDC Short Description	CDC Full Vaccine Name
	198	DTP-hepB-Hib Pentavalent Non-US	Diphtheria, pertussis, tetanus, hepatitis B, Haemophilus Influenza Type b, (Pentavalent)
	200	influenza, Southern Hemisphere, pediatric, preservative free	influenza, seasonal, Southern Hemisphere, quadrivalent, pediatric 0.25mL dose, preservative free
	201	influenza, Southern Hemisphere, preservative free	influenza, seasonal, Southern Hemisphere, quadrivalent, 0.5mL dose, no preservative
	202	influenza, Southern Hemisphere, quadrivalent, with preservative	influenza, seasonal, Southern Hemisphere, quadrivalent, 0.5mL dose, with preservative

## Bug Fix Summary

	This release does not contain bug fixes.

## Product Documentation

Product documentation is located on the STC Documentation Portal:  
<https://documentation.stchome.com/>.

The following documents are available for this version of ImmuCast:

- ImmuCast 5.18.8 User Guide
- ImmuCast 5.23.0 Release Notes