



PHC Hub Release Notes

v5.23.0



Support Services

For general support on this product, contact your system administrator or help desk. For up-to-date documentation, visit the STC Documentation Portal at <https://documentation.stchome.com/>.

Connect with Us on Social Media



Copyrights and Trademarks

© 2019 by Scientific Technologies Corporation (STC). All rights reserved.

This documentation describes the following: PHC Hub (v5.23.0) release notes

No part of this publication may be altered, reproduced, transmitted, stored in a retrieval system, or translated into any human or computer language, in any form or by any means, electronic, mechanical, magnetic, optical, or otherwise, without the prior written permission of the copyright owner.

This document, along with any embedded media, is provided for informational purposes only. Scientific Technologies Corporation (STC) provides this material "as is" and does not assert that this document is error free. The information contained herein may not necessarily accurately represent the current or planned functions of the application, and may be subject to significant and frequent modification. Scientific Technologies Corporation (STC) may change this document, any embedded content, and the product described herein at any time. Any changes will be incorporated in new versions of this document.

Scientific Technologies Corporation (STC) recognizes the rights of the holders of all trademarks used in its publications.

This document may provide hyperlinks to third-party websites or access to third-party content. Links and access to third-party sites are provided for your convenience only. Scientific Technologies Corporation does not control, endorse, or guarantee third-party content and is not responsible for any content, associated links, resources or services associated with a third-party site. Scientific Technologies Corporation shall also not be liable for any loss or damage associated with your use of any third-party content.

Table of Contents

Support Services	2
Connect with Us on Social Media	2
Copyrights and Trademarks	2
Table of Contents	3
New Features/Improvements	4
Known Issues.....	5
Export to Excel.....	5
Product Documentation.....	5

New Features/Improvements

The following are new features, functionalities or enhancements in this version of PHC Hub.

Please note: If you choose to skip this database patch, PHC Hub will still superficially appear to operate normally. However, two consequences will occur:

1. Missing patient data resulting from the lack of the newly added nationality column could affect reports that query cp_hold_patient.
2. The issue resolution default will actually change from Warn to Ignore, changing the behavior of existing import profiles.

Key	Summary	Description
CONSORT-4959 CONSORT-4961	Addition of new Import Profile Issue Resolution for Patient Nationality	The following new issues have been added to PHC Hub Import Profile Issue Resolution Patient Nationality ...is missing (CONSORT-4959) ...is invalid (This change does not require client testing as only inactive codes will trigger this issue and all country codes all active). ...is unrecognized (CONSORT-4961)
CONSORT-4962 CONSORT-4963 CONSORT-4964	Addition of new Import Profile Issue Resolution for Patient Birth Country	The following new issues have been added to PHC Hub Import Profile Issue Resolution Patient Birth Country ...is missing (CONSORT-4962) ...is invalid (included for consistency) (CONSORT-4963) ...is unrecognized (CONSORT-4964)
CONSORT-4958	Newly added issue resolutions are set to ignore by default.	
CONSORT-4965	Nationality and Birth Country Message Logging	As a result of adding these issue resolutions, Nationality and Birth Country messages are logged in the Import Log, Import Profiles > Profile and Batches.
CONSORT-4966	Nationality in queries	PID-28 Nationality is now accepted in updates and output in queries.

Key	Summary	Description
		<p>For example:</p> <pre>PID 1 10069^^^^SR~INT 20.00 ROCK^YUSEF^^^^L 20110611 M 2106- 3^White^CDCPHINVS 20 INTERNATIONAL CIR^^PHOENIX^ARIZONA^85020^USA^M^^MARI COPA 2135-2^Hispanic or Latino^HL70189 N 1 USA^United States^HL70399 N </pre>

Known Issues

The following are known issues with this release:

- The Attestation Report is currently not working.
- Version numbers do not display at the bottom of the left navigation menu.
- As previously noted in the PHC Hub v. March 2018 Release Notes, code changes are required in IWeb to complete the downstream functionality to support the relaxed consistency logic for matching VFC Eligibility and Funding Source. Clients should see fewer message rejections related to PHC Hub processing for VFC Eligibility and Funding Source, but may not see a change in lot decrementing or IIS data retention until changes have been made in IWeb.

Export to Excel

Follow these steps to export the testing steps to an Excel spreadsheet:

1. Click this link and log in to Jira if required:
<https://stchome.atlassian.net/issues/?filter=52713>.
2. Click the **Change View** icon and select **List View**.
3. Click the **Export** icon (it looks like a download icon) and export as needed.

Product Documentation

Product documentation is located on the STC Documentation Portal:
<https://documentation.stchome.com/>.

The following documents are still applicable for this version of PHC Hub:

- PHC Hub (v. July 2018) User Guide (no changes for this version)
- PHC Hub (v. March 2018) Quick Reference Guides (QRGs) – (no changes for this version)