



# VOMS

## Known Issues

V2.16.12.1.1



# Table of Contents

- Define Drop-Down lists search bar not working properly .....2
- Unable to save Cold Storage unit information when existing 'Calibration Expiration' date is removed .....3
- Cannot upload Temp Log File .....4
- On Properties page clicking 'Cancel' when trying to navigate away from page prior to saving does not re-load page .....5
- Filling in then clearing min/max temperature still tries to validate as 0 .....6
- VTrckS Provider Data Files .....7
- Column sorting is dysfunctional on Transfer page when entering in quantity .....8
- Returns are being displayed on the Transfers for Search History .....9
- Min/Max Cold Storage isn't working correctly. .... 10
- On Define Drop-Downs List the 'Presentation' list is showing incorrect information ..... 11

# Define Drop-Down lists search bar not working properly

**Description:** When typing in the search bar on the Define Drop-Down Lists page the page is not filtering based upon the entered text.

**Version Affected:** VOMS v2.16.12.1.1

**Who Identified the Issue:** STC

**Priority:** Medium

**Date:** 03/14/2017

## **Steps to Test:**

1. Login into VOMS
2. Go to Admin & Settings
3. Navigate to Define Drop-Down Lists
4. Start typing in the Search Bar
5. Notice the page is not filtering based on text entered

**Expected Results:** The page will filter based on the text entered in the Search Bar.

**Actual Results:** The page is not filtering based on the text entered in the Search Bar.

# Unable to save Cold Storage unit information when existing 'Calibration Expiration' date is removed

**Description:** When editing information for a Cold Storage Unit the page should successfully save upon clicking Save if the Calibration Expiration is removed.

**Version Affected:** VOMS v2.16.12.1.1

**Who Identified the Issue:** STC

**Priority:** Medium

**Date:** 03/13/2017

## Steps to Test:

1. Login into VOMS
2. Select a Facility
3. Click on Cold Storage menu
4. Click the 'Edit' icon by a Cold Storage unit
5. Un-select the Date of Last Calibration
6. Select Save
7. Notice the Warning message displayed

**Expected Results:** When deselecting the Date of Last Calibration, the page should save successfully IF all required fields are completed.

**Actual Results:** When deselecting the Date of Last Calibration, a warning message will appear upon clicking Save.

# Cannot upload Temp Log File

**Description:** The Import File functionality is not allowing users to import files once all page fields have been cleared.

**Version Affected:** VOMS v2.16.12.1.1

**Who Identified the Issue:** STC

**Priority:** Medium

**Date:** 03/14/2017

## Steps to Test:

1. Login into VOMS
2. Select a Facility
3. Navigate to Cold Storage menu
4. Select Upload Temp Log
5. Select Cold Storage unit
6. Select a Defined Profile
7. Select File and click Import
8. Click Clear button
9. Re-select Cold Storage unit and select Defined Profile
10. Try to Import the same file again
11. System fails in allowing you to Import File
12. Notice you cannot Import a File once you've cleared the fields.

**Expected Results:** A file should Import each time a file is selected and Import is clicked.

**Actual Results:** Once a file is selected, it doesn't populate the file name in the Choose File box after selecting a file from your computer.

**Note:** Two workarounds were identified:

- Select a different file then select the original file. Notice the file is selected
- Navigate away from the page, return to the Upload temperature page and continue with the Import process

# On Properties page clicking 'Cancel' when trying to navigate away from page prior to saving does not re-load page

**Description:** In the Admin & Settings menu when a user is on the Properties page, changes any field and navigates away before saving and clicks 'Cancel' on the warning message the page does not automatically refresh.

**Version Affected:** VOMS v2.16.12.1.1

**Who Identified the Issue:** STC

**Priority:** Medium

**Date:** 03/14/2017

## **Steps to Test:**

1. Login into VOMS
2. Click Admin & Settings
3. Navigate to Properties
4. Change any field
5. Click Home menu
6. Warning messages displays, click Cancel
7. Notice the page does not re-load

**Expected Results:** Upon clicking Cancel in the warning message pop-up, the page should re-load

**Actual Results:** The page does not re-load once Cancel is clicked in the warning message pop-up

# Filling in then clearing min/max temperature still tries to validate as 0

**Description:** When entering temperatures and then clearing the field, users will still see a warning message.

**Version Affected:** VOMS v2.16.12.1.1

**Who Identified the Issue:** STC

**Priority:** Medium

**Date:** 03/15/2017

## **Steps to Test:**

1. Login into VOMS
2. Click on Cold Storage
3. Switch to Min/Max
4. Enter a number into the Min field
5. Clear the entry
6. Enter a value below '0' for the Max field
7. Notice a warning message in the Min field

**Expected Results:** No warning message will display if entry is cleared.

**Actual Results:** A warning message will display and the user will not be able to save without entering a value.

**Note:** This can be worked around if you enter in a value that is less than the Max value entered.

# VTrckS Provider Data Files

**Description:** Sometimes the search parameters and downloading of the Provider Data Files are failing.

**Version Affected:** VOMS v2.16.12.1.1

**Who Identified the Issue:** STC

**Priority:** Medium

**Date:** 03/15/2017

## **Steps to Test:**

1. Login into VOMS
2. Go to Admin & Settings
3. Click on VTrckS Provider Master Data Export
4. Select organization that has an approved Provider Agreement
5. Select check box for "Only include orgs/facilities with provider agreement in the 'Approved' Status"
6. Notice that there are no results when there should with one or more providers with an approved Provider Agreement

**Expected Results:** Each time a parameter is defined the page should load accordingly and the file should successfully download.

**Actual Results:** On occasion the search parameters are not filtering the page correctly and the Provider Files are failing to download.



# Column sorting is dysfunctional on Transfer page when entering in quantity

**Description:** Table header sort function is not working as intended upon entering in quantity to transfer.

**Version Affected:** VOMS v2.16.12.1.1

**Who Identified the Issue:** STC

**Priority:** Medium

**Date:** 03/16/17

## **Steps to Test:**

1. Login into VOMS
2. Navigate to Orders & Returns and select Orders & Transfers
3. Click New transfer
4. Sort a column (Any column except vaccine)
5. Enter in quantity
6. System automatically removes the sort and you are shown that entered in quantity for the first available vaccine but cannot comment on this amount.
7. Sort again with the same column initially chosen and your quantity comes back to the correct vaccine and you can then comment.

**Expected Results:** The table should sort and save the original entered quantity

**Actual Results:** The table sorting functionality changes once a quantity is entered.

**Note:** A work around is to not use the table header sorting functionality.

# Returns are being displayed on the Transfers for Search History

**Description:** On the Search History table under the Transfers tab, Returns are also being shown.

**Version Affected:** VOMS v2.16.12.1.1

**Who Identified the Issue:** STC

**Priority:** Medium

**Date:** 03/17/2017

## Steps to Test:

1. Login into VOMS
2. Login to the VOMS
3. Select an Org/Fac
4. Navigate to the Search History screen.
5. Select Transfers tab
6. Click Search.
7. Notice Returns are displayed along with Transfers.

**Expected Results:** Only Transfers should be displayed on the Transfer tab

**Actual Results:** Both Transfers and Returns are being displayed on the Transfer tab

# Min/Max Cold Storage isn't working correctly.

**Description:** Cold Storage Export files are incomplete.

**Version Affected:** VOMS v2.16.12.1.1

**Who Identified the Issue:** STC

**Priority:** Medium

**Date:** 03/16/2017

## Steps to Test:

1. Login into VOMS
2. Go to Cold Storage and Upload a Temperature Log
3. Submit Temp log to State as 2X Daily Temps
4. Select as of date within the file date
5. Export data
6. Notice that the file it shows the correct date range but it does not contain any data from the temp log file.

**Expected Results:** The exported Temperature File should contain all dates and temperature results uploaded.

**Actual Results:** The exported Temperature File may not contain all Temperature data and should not be relied on.

# On Define Drop-Downs List the 'Presentation' list is showing incorrect information

**Description:** On Define Drop-Down lists, the 'Presentation of the vaccine vials, syringes etc' is displaying funding source information instead of the presentation information.

**Version Affected:** VOMS v2.16.12.1.1

**Who Identified the Issue:** STC

**Priority:** Medium

**Date:** 03/17/2017

**Steps to Test:**

1. Login into VOMS
2. Navigate to Admin & Settings > Define Drop-Downs List.
3. Select Presentation of the vaccine vials, syringes etc.
4. Notice the Funding source values are displayed.

**Expected Results:** Presentation information should be displayed.

**Actual Results:** Funding source values are being displayed.